

- ✓ Mid-Turbinate Nasal Swab
- ✓ Easy to Follow Collection Instructions
- ✓ Online Registration
- ✓ Quick Online Results
- ✓ FDA EUA Approved

About Our Test

This test detects the presence or absence of SARS-CoV-2, the virus that causes COVID-19¹. Sample collection is an easy process using a nasopharyngeal swab. Please note this test does not detect antibodies or immunity; it is designed to determine active infection with the virus that causes COVID-19.

How It Works

- Register the test packet online.
- Collect and return the sample to our lab.
- Results typically available on our secure patient portal 24-48 hours after receipt to our lab.


Why Our Test?

The “mid-turbinate” swab is less invasive than the deep nasal pharyngeal and easier to administer.

This test is run in our CLIA (Clinical Laboratory Improvement Amendments act 1988.) certified and CAP (College of American Pathologists) accredited High Complexity Laboratory. Sensitivity and specificity is the 98-99%. This test is FDA approved under Emergency Use Authorization and is CLIA certified.

What is Included?

The SARS-CoV-2 Test Collection packet includes: Requisition form, symptoms questionnaire, collection and results instructions, NP swab and collection tube, bar codes, return zip lock bag and FedEx return clinical pak.



Item # 1005
MSRP: \$129.99

¹Our lab will test for the presence of three SARS-CoV-2-specific genes by reverse transcriptase real-time PCR (rRT-PCR; TaqPath COVID-19 Combo Kit, Thermo Fisher). The NP Screening Test for SARS-CoV-2 uses a sensitive, real-time RNA-PCR assay specifically designed for the detection of severe acute respiratory syndrome Corona virus 2 (SARS-CoV-2).